



PROTECTING WHAT MATTERS

# Pig DOT Packaging Instructions


1



**ERECT BOX**

Pull out sides and push bottom down until tab engages


2



**PUT IN BOTTOM FOAM**

Place bottom foam into the base of the box


3



**LINERS & FOAM DISC**

Ensure liner and foam disc are in each pig base


4



**INSERT VIALS**

Place product filled vial into pig base and put on lid

5



**TAPE LID**

Tape lid to base all the way around


6



**SEAL INTO BAG**

Put pig into plastic bag and seal


7



**SET PIGS INTO FOAM**

Put desired number of bagged pigs into bottom foam


8



**ADD ICE PACKS (OPTIONAL)**

Add ice packs to the 4 outer slots in the lower foam base if needed


9



**PUT IN TOP FOAM**

Add top foam piece securely to cover the pigs

10



**PLACE TRACKER (OPTIONAL)**

Add tracking device to top insert of foam cover

11



**CLOSE BOX**

Close top of box and engage flap and locking tab



Please contact us if you have any questions.